

Firework™ Stiff Guidewires

Product type	Model & Specifications	Tip form	Rod form	Diameter	Length	Crimp diameter of head	Packaging
Firework - Amplatz Extra-stiff guidewire	SSG-35-2750-S	Straight	Wire coil cladding	0.035 in/0.89mm	2750mm		5 Pkg/Box
	SSG-35-2750-J	J-type				3mm	5 Pkg/Box
Firework - Safari Super-stiff guidewire	SSG-35-2750-SC	Double curve (small)				32mm*29mm	5 Pkg/Box
	SSG-35-2750-LC	Double curve (large)				42mm*42mm	5 Pkg/Box
Firework - Lunderquist Super-stiff guidewire	SSG-35-2750-ES	Straight	Polished rod body				5 Pkg/Box
Firework -Super Stiff Safari Double-curve guidewir	SSG-35-2750-ESC	Double curve (small)				32mm*29mm	5 Pkg/Box
	SSG-35-2750-ELC	Double curve (large)		42mm*42mm	5 Pkg/Box		