

Magpie™ Balloon Dilatation Catheter

| Model | Balloon diameter Ø (mm) | Outer diameter of catheter | Minimum sheath dimensions | Nominal pressure (atm) | Rated burst pressure (atm) |
|---------|-------------------------|----------------------------|---------------------------|------------------------|----------------------------|
| VB1240L | 12 | 7F | 8F | 4 | 7 |
| VB1440L | 14 | 7F | 8F | 4 | 7 |
| VB1640L | 16 | 8F | 12F | 4 | 7 |
| VB1840L | 18 | 8F | 12F | 4 | 7 |
| VB2040L | 20 | 9F | 12F | 2 | 6 |
| VB2240L | 22 | 9F | 14F | 2 | 5 |
| VB2340L | 23 | 9F | 14F | 2 | 5 |
| VB2540L | 25 | 9F | 16F | 2 | 4 |
| VB2640L | 26 | 9F | 16F | 2 | 4 |
| VB2840L | 28 | 11F | 16F | 2 | 4 |