

Stiff Guidewires

Product type	Model & Specifications	Diameter	Length	Head	Recommended use
Vigo Track™ -stiff guidewire	SSG-35-2750-ES	0.035 in	2750mm	Straight head	TEVAR/EVAR
Safari Super-stiff guidewire	SSG-35-2750-SC			Small (32mm*29mm)	TAVR
	SSG-35-2750-LC			Large (42mm*42mm)	
Super Stiff Safari Double-curve guidewir	SSG-35-2750-ESC			Small (32mm*29mm)	
	SSG-35-2750-ELC			Large (42mm*42mm)	