

Stiff Guidewires					
Brand name	Specification	Diameter	Length	Head	Recommended use
VigoTrack™	SSG-35-2750-ES	0.035 in	2750mm	Straight	TEVAR/EVAR
Archmid™	SSG-35-2750-SC			Small (32mm*29mm)	TAVR
	SSG-35-2750-LC			Large (42mm*42mm)	
Archmid-E™	SSG-35-2750-ESC			Small (32mm*29mm)	
	SSG-35-2750-ELC			Large (42mm*42mm)	