

Magpie™ Valvuloplasty Balloon Dilatation Catheter

Specification	Balloon diameter(mm)	Catheter outer diameter	Recommended sheath size	Nominal pressur(atm)	Rated burst pressure(atm)
VB1240L	12	7F	8F	4	7
VB1440L	14	7F	8F	4	7
VB1640L	16	8F	10F	4	7
VB1840L	18	8F	10F	4	7
VB2040L	20	9F	12F	2	6
VB2240L	22	9F	12F	2	5
VB2340L	23	9F	14F	2	5
VB2540L	25	9F	14F	2	4
VB2640L	26	9F	16F	2	4
VB2840L	28	11F	16F	2	4