

DragonFly Transcatheter Mitral Valve Repair System			
Components Name	Specification	Sheath O.D.	Sheath length (mm)
Guiding Sheath	GS24800	24	800
Components Name	Specification	Clip arm width(mm)	Clip arm length (mm)
Valve Clip System	SN0409	4	9
	XN0412	4	12
	SW0609	6	9
	XW0612	6	12
Reusable accessories are intended to be used with DragonFly System product family during the DragonFly procedure. DragonFly Accessories: Stabilizer(DSS-45-45-20-S、DSS-45-45-20L)、Table(DSS-T)、Support Plate(DSS-SP)			